

MODEL DH4-38GD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	38 Cu. Ft.
Capacity [Liters]	1060 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F
Dimensions	Width x Depth x Height
Exterior [Inches]	47.25" x 29.25" x 78.5"
Interior [Inches]	43" x 23.25" x 60.5"
Exterior [Centimeters]	120 x 74 x 199
Interior [Centimeters]	109 x 59 x 153
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	2 Doors, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter (Top Back Wall)
Casters	4, Locking
Storage	10 Adjustable Shelves
Shipping Weight	500 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 8.3A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts

Warranty

3 Years Parts & Labor.

Additional 4 Years on Compressor Parts, No Labor.



IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

