

MODEL DHH20-14MDP

MANUAL DEFROST FREEZER -15°C TO -25°C

Capacity [Cu. Ft.] 14 Cu. Ft.
Capacity [Liters] 396 Liters
Temperature Range [°C] -15°C to -25°C
Temperature Range [°F] 5°F to -13°F

Dimensions Width x Depth x Height
Exterior [Inches] 23.63" x 24.38" x 73"
Interior [Inches] 19" x 17.5" x 63.75"
Exterior [Centimeters] 60 x 62 x 185
Interior [Centimeters] 48 x 44 x 161

Compressor Hermetic Compressor
Air Circulation Gravity Flow
Defrost Method Manual
Refrigerant R290
Exterior Powder Coating
Interior Molded ABS
Insulation Polyurethane Foam
Lockable Door Keyed
Porthole 3/4" Diameter
Leveling Legs 4
Storage 6 Fixed Shelves
Shipping Weight 225 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 3.5A
NEMA Configuration NEMA 5-15
Certification(s) ETL Listed
Energy Star Rated

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp
Alarm Relay Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Year Parts & Labor