



Lightweight and portable, the MTB baby and toddler scale offers excellent value and versatility for pediatricians and other healthcare providers. With a capacity of 20kg, this scale boasts the flexibility to weigh newborn babies to toddlers on the same device. MTB features quick start-up and easy operation, so minimal training is needed. A cleverly designed ergonomic cradle locks securely in place and provides comfortable support for babies. The cradle removes easily to reveal a non-slip rubber mat underneath, simplifying the task of weighing toddlers. A memory function retains the most recent weight, allowing the user to record results after safely removing the child from the scale.

Offering the flexibility to weigh newborns and toddlers on the same scale



Features

- High-contrast display is easily visible from a distance
- Memory feature can recall the last weight taken, allowing the user to focus on the child
- Cradle can be removed to weigh toddlers on the non-slip rubber mat
- Ergonomic cradle provides comfortable support
- Lightweight design simplifies transport and portability
- External calibration allows for verification and adjustment with weights
- Simple operation allows quick startup and use with little or no training
- The auto power-off feature conserves battery life
- Operates with AC adapter (included) or batteries (not included)

Applications

- Weighing

About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

MTB Baby and Toddler Scale

Model	MTB 20
Capacity	44lb / 20kg
Readability	0.005lb / 5g
Repeatability (S.D.)	0.01lb / 5g
Linearity (+/-)	0.01lb / 10g
Pan Size	11.2"x22" / 285x560mm
Weighing Units	kg, lb
Stabilisation Time (sec)	2
Calibration	External calibration
Display	LCD with 1.1" / 29mm -high digits
Power Supply	6VDC 500mA adapter / 4 x AA batteries
Operating Temperature	10° to 40°C
Housing	ABS plastic
Overall Dim.	22"x13.8"x5.5" / 560x350x140mm (wxdxh)
Net Weight	7.1lb / 3.2kg

Accessories

Item number	Description
700660290	Calibration certificate

Supplied by