MDW Mechanical Physician Scales



Adam Equipment's MDW mechanical scales offer an affordable solution for checking the weight of fitness center or gym members, and patients in hospitals and doctors' offices. Wheels allow easy mobility and a notched counter weight clearly shows the numbers. MDW scales feature rugged construction with a metal base and a non-slip weighing platform to ensure stability and comfort. Most models include an adjustable height rod to simultaneously measure height and weight. The MDW 200B does not include wheels or height rod.

HEADAM

Affordable, easy-to-operate scale measures height and weight



Features

- Height measurements from 60cm to 212cm
- Integrated wheels permit easy movement
- · Zero adjustment
- Anti-slip platform ensures safety when stepping on and off the scale
- · Height rod facilitates measuring children and adults while weighing
- · Metal construction for base and pillar

ApplicationsWeighing

About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

MDW Mechanical Physician Scales



Model	MDW 160M	MDW 200B	MDW 200M
Capacity	350lb / 160kg	440lb / 200kg	440lb / 200kg
Readability	0.2lb / 0.1kg		
Repeatability (S.D.)	0.4lb / 0.2kg		
Linearity (+/-)	0.4lb / 0.2kg		
Pan Size	14.8"x10.8" / 375x275mm		
Weighing Units	kg, lb		
Stabilization Time (sec)	3		
Height Rod	Yes	No	Yes
Calibration	N/A		
Display	Graduated beam		
Operating Temperature	5° to 40°C / 41° to 104°F		
Overall Dim.	20.9"x10.8"x58.5" / 530x275x1485mm (wxdxh)		
Net Weight	33.1lb / 15kg		

Supplied by

www.adamequipment.com

Adam Equipment follows a policy of continuous improvement and reserves the right to change specification and pricing without notice. Adam Equipment's standard terms and conditions of sale apply. Errors and omissions excepted.

 © 2017 Adam Equipment
 EN-6062
 03/23/17