



Part Number	Description	Code	Remarks
281395	Protective cover for MFD	OA154	Silicon protective cover protects the main unit from dirt and scattering of samples
281396	Protective cover for MFH	OA155	
281394	External temperature sensor	OA153	Sensor used when external temperature control is performed with MFH
281381	Pole set	OA143	Pole for fixing baths and mixers $\Phi 10 \times 480$ mm
281384	Container fall prevention frame	OA146	Variable frame to prevent containers such as beakers from falling when placed on the stirring table
281385	Stage for lab jack	OA147	Allows main unit to be placed on a lab jack
281587	Power cord (round terminal 2m)	OA183	Round terminal power cable
281382	Aluminum block handle	OA144	Handle to carry hot aluminum block with one hand
281383	Dual-handed aluminum block handle	OA145	Handles to carry hot aluminum blocks with both hands
231632	Muff	OLM44	$\Phi 5 \sim \Phi 13$ mm
231633	Muff	OLM46	$\Phi 6 \sim \Phi 17$ mm
231634	Muff	OLM48	$\Phi 9.5 \sim \Phi 29$ mm
231635	Double opening clamp	OLM50	Tightening adjustment range 3~55 mm, shaft $\Phi 10$ mm
231636	Double opening clamp	OLM52	Tightening adjustment range 3~80 mm, shaft $\Phi 20$ mm
222193	Glass tank for BOG100	OBO14	Mounted on MFH and used as an oil bath $\Phi 150$ mm 1.0L
222194	Glass tank for BOG200	OBO16	Mounted on MFH and used as an oil bath $\Phi 180$ mm 2.2L
281386	High magnetic agitator	OA148	Oval $\Phi 6 \times 15$
281390	High magnetic agitator	OA149	Octagon $\Phi 3 \times 13$
281391	High magnetic agitator	OA150	Octagon $\Phi 8 \times 13$
281392	High magnetic agitator	OA151	Octagon $\Phi 8 \times 38$
281393	Magnetic agitator	OA152	Micro $\Phi 2 \times 5$, set of 5
F-4028-02	Magnetic agitator	TB-20	$\Phi 7 \times 20$ 12 pcs.
F-4028-03	Magnetic agitator	TB-30	$\Phi 8 \times 30$ 12 pcs.
F-4028-04	Magnetic agitator	TB-40	$\Phi 8 \times 40$ 12 pcs.
F-4025-04	Magnetic agitator	A-43	$\Phi 13 \times 43$ 6 pcs.




Protective cover (MFH)

281396




Stage for lab jack

281385




Container fall prevention frame

281384




Dual-handed aluminum block handle

281383



Aluminum block handle

281382



High magnetic agitator

281390 / 281391 / 281392