

SINCE 1889



STANDARD LABORATORY USE STERILIZER

SN200C/210C/300C/310C/500C/510C

The SN Sterilizers are ergonomically designed and feature easy top loading and handling of samples

Features

- Comfortable accessibility and easy operation
- Easy to clean and maintain
- Mobile on wheels
- With high performance and safety functions
- Programmable or preset programs for commonly used sample types
- Reliable high pressure steam sterilization up to 135°C
- Cooling fan to shorten cool down time
- Various add-on options and wide selection of baskets/containers available



Specifications

Model	SN200C/210C	SN300C/310C	SN500C/510C
Operating temp. range	45~135°C		
Max. operational pressure	0.255 MPa		
Cooling	Equipped with cooling fan		
Safety devices	Self diagnostic function, locking lid mechanism, safety pressure valve		
Internal dimensions	Dia 300 x 305 mm	Dia 300 x 445 mm	Dia 300 x 665 mm
Internal capacity	20L	32L	47L
Power source AC100-120V 50/60Hz Power source AC200-240V 50/60Hz	SN200C 20~16.5A SN210C 10~8.5A	SN300C 20~16.5A SN310C 10~8.5A	SN500C 23.5~19.5A SN510C 12~10A
Weight	~65kg	~75kg	~85kg