# **Semi-automatic Benchtop Glassware Washer**

## AW47-115V / AW47-220V



### Washing time Setting range 0~60 min.

# Easy to use benchtop semi-automatic glassware washer

- Semi-automatic washer, easy to operate by simply setting time then start
- Upward and downward two-way pressurized water jet method with rotating jet nozzles bring high level cleaning. Detergent washing is also available
- With built-in water heater, no boiler piping and water heating system are required
- Optional jet rack is available for hard to clean targets, such as glassware with narrow neck or body

### **Control Panel**



### Specifications

#### AW47-115V Model AW47-220V Cleaning method Upward and downward two-way pressurized water jet method Rotating jet nozzles (fixed when using jet rack) Washing water temp. Room temp. ~ 60°C Built-in heater 1kW, room temperature to 60°C Water heater 0.1~0.3MPa Supply water pressure Glassware stand Table (Standard), racks (optional) Water supply Electromagnetic valve open/close Water drain Natural drainage by water level gap Exterior material Chrome-free electric galvanized steel plate, chemical-resistant paint Interior material Stainless steel W450 x D490 x H875mm External dimensions Internal dimensions W420 x D450 x H570mm Pump 200W Spin table Dia. 420mm Drop down style Door Weight ~43kg AC115V 13A Power source (50/60Hz) AC220V 7A Included accessories Water supply hose (with coupler) 2m 1pc. Drain hose (I.D.25.4mm) 1.5m 1pc. Vinyl cover 1pc. Phosphorus-free detergent 1kg (50ml measuring spoon 1pc.) Water supply unit 1set Consumable Phosphorus-free detergent

#### \* Protrusions excluded







Jet rack (glassware not included)



Test tube rack (glassware not included)



Phosphorus-free detergent

No.	Product name	Description	Product code
(1)	Jet rack	Hold up to 36 pcs. of 100ml flask	291090
(2)	Test tube rack	Hold up to 450 pcs. of ø18.5mm test tube	291091
(3)	Detergent	Phosphorus-free detergent 8kg	8190026001