



*96-LED ring polarized light illuminator,
multi-plug*

Chromaticity and LED	Coordinates (x,y)	(0.32, 0.32)
	Color Temperature (K)	6,300
	LED number	96
	Brightness control	Manual
	Lifetime (hours)	> 65,000
	Max. required power (W)	6
Optical power	Luminous flux (lumen)	480
	Illuminance (at 10 cm distance, lux)	9,000 (on 6 cm diameter)
Power Supply for Illumination	Type	External
	Connector	Jack, 1.4 mm
	Power plug type	Multi-plug (EU, UK, US)
	Input voltage	100/240 Vac, 50/60 Hz
	Output voltage	24 Vdc 1.5 A
Accessories Included	Fixing ring	Yes
	Rotating Polarizing filter	Yes
	Analyzer filter	Yes