



*Coaxial coarse & fine focusing system*

<b>Focusing System</b>	Type	Coaxial coarse & fine
	Coarse total travel (mm)	50
	Fine total travel (per single rotation) (mm)	0.2
	Fine graduations	100
	Fine resolution ( $\mu\text{m}$ )	2
	Adjustable tension	Yes
	Head holder internal diameter (mm)	76
<b>Accessories Included</b>	User Manual	Digital version (downloadable)
<b>Product Dimensions</b>	Total height (mm)	100
	Total width (mm)	185
	Total depth (mm)	250
<b>Product Weight</b>	(kg)	1.1