

# TW 3XTL

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	3XTL
LN <sub>2</sub> Capacity (L)	3
Static Hold Time (Days)	27
Evap. Rate (Liters/Day)	0.11
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	7.6
Height (in)	17.2
WEIGHT	
Empty (Lbs)	7.2
Full Weight (Lbs)	12.5
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	N/A



800.886.2796 | Fax 317. 297.7988  
sales@ICBiomedical.com  
www.ICBiomedical.com