

TW 12HCL

LN₂ HIGH-CAPACITY SMALL FREEZER



For long-term storage of large quantities at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

| SPECIFICATIONS | 12HCL |
|--|-------|
| LN ₂ Capacity (L) | 12 |
| Static Hold Time (Days) | 60 |
| Evap. Rate (Liters/Day) | 0.2 |
| DIMENSIONS | |
| Necktube (in) | 3.6 |
| Overall Dia. (in) | 15.6 |
| Height (in) | 19 |
| WEIGHT | |
| Empty (Lbs) | 21.6 |
| Full Weight (Lbs) | 43 |
| CAPACITY | |
| Number of Canisters | 6 |
| Number of .5 cc Straws – Bulk (1 Level) | 2940 |
| Number of .25 cc Straws – Bulk (1 Level) | 6060 |
| Number of 1.2 & 2.0 ml Vials | N/A |



800.886.2796 | Fax 317. 297.7988
 sales@ICBiomedical.com
 www.ICBiomedical.com