

# TW 20XT

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER



For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	20XT
LN <sub>2</sub> Capacity (L)	20.7
Static Hold Time (Days)	230
Evap. Rate (Liters/Day)	0.09
DIMENSIONS	
Necktube (in)	2.00
Overall Dia. (in)	15.6
Height (in)	25.8
WEIGHT	
Empty (Lbs)	26
Full Weight (Lbs)	62.8
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	150



800.886.2796 | Fax 317. 297.7988  
[sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)