

# TW LS6000

## VIAL STORAGE



For liquid- or vapor-stage vial storage in convenient box-type racks.



SPECIFICATIONS	LS6000
LN <sub>2</sub> Capacity (L)	165
Static Hold Time (Days)	194
Evap. Rate (Liters/Day)	0.84
DIMENSIONS	
Necktube (in)	8.5
Overall Dia. (in)	26.9
Height (in)	39
WEIGHT	
Empty (Lbs)	121
Full (Lbs)	410

## FEATURES

- Convenient Storage Systems – rack index location ring and internal location fixture, computer compatible box storage for simple inventory management
  - Number of Racks: 6
  - Number of Shelves / Rack: 10
  - 2ml Vial Capacity: 6000
- In vapor-phase storage, temperature at the top of racks is -190°C.
- High strength aluminum body and magneformed necktube design
- Superior vacuum performance with super insulation provides maximum holding times
- Security – Accessory low-level alarm and CS100 controller are available with remote connection capability
- Lockable lids



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)

[www.ICBiomedical.com](http://www.ICBiomedical.com)

Made In USA