

TW 25LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	25LD
LN ₂ Capacity (L)	25
Static Hold Time (Days)	109
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in)	2.5
Overall Dia. (in)	15.6
Height (in)	25.8
WEIGHT	
Empty (Lbs)	23.2
Full (Lbs)	67.7

FEATURES

- 25 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA