

TW CLASSIC-25

LN₂ DEWAR



For safe storage and dispensing of small accounts of liquid nitrogen.



SPECIFICATIONS	CLASSIC-25
LN ₂ Capacity (L)	25
Static Hold Time (Days)	119
Evap. Rate (Liters/Day)	0.21
DIMENSIONS	
Necktube (in)	2.0
Overall Dia. (in)	15.5
Height (in)	22.9
WEIGHT	
Empty (Lbs)	19
Full (Lbs)	63.5

FEATURES

- 25 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com

www.ICBiomedical.com

Made In USA

© 2021, IC Biomedical, LLC