

TW 50LD

LN₂ DEWAR



For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	50LD
LN ₂ Capacity (L)	50
Static Hold Time (Days)	122
Evap. Rate (Liters/Day)	0.41
DIMENSIONS	
Necktube (in)	2.5
Overall Dia. (in)	18.8
Height (in)	32.4
WEIGHT	
Empty (Lbs)	38.7
Full (Lbs)	127.7

FEATURES

- 50 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA