

ARIES

VP CARTRIDGES (STANDARD)

LAB SERIES, VAPONICS / BARNSTEAD COMPATIBLE

LAB REPLACEMENT CARTRIDGES

VP Series cartridges provide ultra high purity filtration for a variety of applications. Constructed of non-pigmented, virgin materials the spin welded construction provides the longest life and highest quality. Applications include deionization, submicron filtration, organics, chlorine, oxygen and amine removal.



HIGHLIGHTS

- Spin welded construction
- High purity, non-pigmented plastic components
- Individually lot numbered
- Highest capacity cartridges
- Highest quality medias
- Vaponics® and Barnstead® Compatible
- Cartridge kits available
- Custom filters available

SPECIFICATIONS

- Max Pressure of 30 psi (206.84 kPa)
- Service Flow 1 gpm (3.79 lpm)
- Volume of 134.2 cu.in (.078 cu ft)
- Standard end cap
- Body Length: 16.31" (Overall 17.0")
Diameter: 3.38"
Weight: 4.8 Lbs (Approx.)

Revision 1.0

© 2020 ResinTech, Inc.



CAPACITY	DESCRIPTION	BARNSTEAD® EQUIVALENT	PART NUMBER
1360 grains	Color Changing High Purity Mixed Bed	-	VP-17-4030
1590 grains	Color Changing High Capacity Mixed Bed	-	VP-17-4040
2975 grains	Color Changing Cation - Condensate Polishing	-	VP-17-3002
1360 grains	High Purity Mixed Bed	D0809, D5027, D8816	VP-17-4010
1360 grains	Mixed Bed Low TOC	D0820	VP-17-4011
2000 gallons (7,571 liters)	Activated Carbon	D0813	VP-17-1050
3670 grains	Cation (Hydrogen Form)	D0815	VP-17-3001
2580 grains	Anion (Hydroxide Form)	D0760	VP-17-3203
1980 ppm	Oxygen Removal	D0811, D8817	VP-17-3680
2000 gallons (7,571 liters)	Pretreatment - Organics Removal	D8815	VP-17-4090
2,000 gallons / 490 grains	High Purity - Organics	D0832	VP-17-4205
400 grains / 660 ppm	High Purity / Oxygen Removal	D8809	VP-17-4200
300 grains / 525 ppm	High Purity / Oxy Removal / Org. Removal	D8811	VP-17-4301
1,590 grains	High Capacity	D0803	VP-17-4203
3,635 grains	Softening	-	VP-17-3000
1,300 grains	Two Bed Layered	-	VP-17-4202
1,300 grains	Four Bed Layered	-	VP-17-4402
2,000 gallons	Condensate Feed Back (Organic Removal)	-	VP-17-4090-OR
3,200 gallons	Condensate Feed Back (Amine Removal)	-	VP-17-3001-AR
1,000 gallons	Condensate Feed Back (Demineralization)	-	VP-17-4010-DI
8,000 gallons	Carbon with 1 lb KDF	-	VP-17-2001
680 grains	Half Sized - High Purity Mixed Bed	D50313	VP-08-4010
1290 grains	Half Sized - Anion (Hydroxide Form)	D50214	VP-08-3203
1,000 gallons / 245 grains	Half Sized - Pretreatment - Organics	D50215	VP-08-4205
1,000 gallons / 245 grains	Half Sized - High Purity / Organics Removal	D50216	VP-08-4205
200 grains / 330 ppm	Half Sized - High Purity / Oxygen Removal	D50217	VP-08-4200
-	VP Series Kit for RO/DI Feed - (3) High Purity Cartridges	-	VPK-4010
-	VP Series Kit for Tap Water Feed - (1) Organics pretreatment & (2) High Purity Cartridges	-	VPK-3805

Revision 1.0

© 2020 ResinTech, Inc.

