

DATA SHEET E 218

DIAPHRAGM VACUUM PUMP N 920 G



Technical features

- 100 % oil-free transfer thereby pure transferring and evacuation of gases
- High suction speed, particularly in the low vacuum range
- Integrated gas ballast valve
- Integrated rotational speed control enables pumping capacity to be easily adapted manually to process requirements
- Ideal for aggressive/corrosive gases and vapors
- Maintenance-free
- Environmentally friendly

Technical data

	4.00	
Flow rate [m ³ /h] at atm. pressure	1.26	
Flow rate [I/min] at atm. pressure	21	
Ultimate vacuum [mbar abs.]	2.0 < 5.0 with gas ballast valve in use	
Operating pressure [bar]	0.5	
Permissible media temperature [°C]	+5 +40	
Permissible ambient temperature [°C]	+10 +40	
Mains [V/Hz]	100 - 240/50 - 60	
Motor protection	IP 20	
Power P ₁ [W]	135	
Operating current [A]	1.4	
With overcurrent protection and main fuse		

OTHER DATA

Hose connections [mm]	for tube ID 10
Weight [kg]	8.5
Dimensions W x H x D [mm]	158 x 226 x 324

MATERIALS

Pump head	PPS
Diaphragm	PTFE-coated
Valves	FFPM

Dimensions



ACCESSORIES	ID
Silencer G 1/8	007006
Fine control valve with vacuum gauge	112432
Small flange, stainless steel SF 16	046625
Spares kit	111905



Performance characteristics





